

## Anti-Pla2G6 Rabbit Polyclonal Antibody For Wb And Elisa - O60733

Item Number: CM0026382

### Introduction

Rabbit polyclonal antibody specific to PLA2G6, validated for western blotting (WB) and ELISA. Cross-reacts with rat. Ideal for research on phospholipid metabolism, mitochondrial function, and neurodegeneration.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PLA2G6 Rabbit pAb                                                                                                    |
| Remarks/Alias    | GVI; PLA2; INAD1; NBIA2; iPLA2; NBIA2A; NBIA2B; PARK14; PNPLA9; Cal-PLA2; IPLA2-VIA; iPLA2beta; PLA2G6               |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 8398                                                                                                                 |
| Gene ID          | 8398                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 151-250 of human PLA2G6 (NP_003551.2). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Rat                                                                                                                  |
| SWISS            | O60733                                                                                                               |
| Protein Weight   | 90kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |