


Anti-Pla2G7 Polyclonal Antibody For Wb, Elisa - Q13093

Item Number: CM0013269



Introduction

Rabbit polyclonal antibody to human PLA2G7 (PAF-AH, LP-PLA2). Validated for WB and ELISA with cross-reactivity to mouse. Use for studying platelet-activating factor degradation and inflammatory oxidative stress pathways.

[Learn More](#)

Parameter	Value
Product Name	PLA2G7 Rabbit pAb
Remarks/Alias	PAFAD; PAFAH; LP-PLA2; LDL-PLA2; PLA2G7
Species	Human
Gene ID (Human)	7941
Gene ID	7941
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 51-240 of human PLA2G7 (NP_005075.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13093
Protein Weight	50kDa
Shipping	Ice bag