


Anti-Pla2G7 Rabbit Polyclonal Antibody For Wb, Elisa - Q13093

Item Number: CM0029337



Introduction

Rabbit polyclonal antibody targeting human PLA2G7 (LP-PLA2, PAF-AH). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying inflammation, oxidative stress, and lipid metabolism in research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	PLA2G7 Rabbit pAb
Remarks/Alias	PAFAD; PAFAH; LP-PLA2; LDL-PLA2; PLA2G7
Species	Human
Gene ID (Human)	7941
Gene ID	7941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 172-441 of human PLA2G7 (NP_005075.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13093
Protein Weight	50kDa
Shipping	Ice bag