

Anti-PLA2R1 Rabbit Monoclonal Antibody For Wb And Elisa - Q13018

Item Number: CM0007087



Introduction

High-quality PLA2R1 rabbit monoclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for phospholipase A2 receptor studies and kidney disease diagnostics research.

[Learn More](#)

Parameter	Value
Product Name	PLA2R1 Rabbit mAb
Remarks/Alias	PLA2R; PLA2-R; PLA2IR; CLEC13C; PLA2G1R; PLA2R1
Species	Human
Gene ID (Human)	22925
Gene ID	22925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 810-910 of human PLA2R1 (NP_031392.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13018
Protein Weight	169 kDa
Shipping	Ice bag