

Anti-Pla2R1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13018

Item Number: CM0016766



Introduction

Rabbit polyclonal antibody against PLA2R1, validated for WB and ELISA in human, mouse, and rat samples. Ideal for studying secretory phospholipase A2 receptor signaling and clearance pathways.

[Learn More](#)

Parameter	Value
Product Name	PLA2R1 Rabbit pAb
Remarks/Alias	PLA2R; PLA2-R; PLA2IR; CLEC13C; PLA2G1R; PLA2R1
Species	Human
Gene ID (Human)	22925
Gene ID	22925
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 870-1130 of human PLA2R1 (NP_031392.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13018
Protein Weight	169kDa
Shipping	Ice bag