

Anti-Placental Alkaline Phosphatase (Plap) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05187

Item Number: CM0008038



Introduction

High-quality rabbit monoclonal antibody targeting human placental alkaline phosphatase (PLAP). Suitable for WB, IHC-P, and ELISA applications, with cross-reactivity to mouse and rat. Supports placental biology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Placental alkaline phosphatase (PLAP) Rabbit mAb
Remarks/Alias	ALP, IAP, ALPI, PALP, PLAP, PLAP-1, Placental alkaline phosphatase (PLAP)
Species	Human
GeneID (Human)	250
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Placental alkaline phosphatase (PLAP) (P05187).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05187
Protein Weight	58kDa
Shipping	Ice bag