

Anti-Placental Lactogen (Csh1) Monoclonal Antibody For Elisa And If-P - P01243

Item Number: CM0002395

Introduction

Rabbit monoclonal antibody targeting human Placental lactogen (CSH1), suitable for ELISA and IF-P. Designed for placental hormone and reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	Placental lactogen (CSH1) Rabbit mAb
Remarks/Alias	PL; CSA; CS-1; CSMT; GHB3; hCS-1; hCS-A; Placental lactogen (CSH1)
Species	Human
Gene ID (human)	1442
Gene ID	1442
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Placental lactogen (CSH1) (CSH1) (P0DML2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	P01243
Protein Weight	25kDa
Shipping	Ice bag