


Anti-Placental Lactogen (Csh1) Polyclonal Antibody For Ihc-P, Elisa - P0Dml2

Item Number: CM0022401



Introduction

Placental lactogen (CSH1) rabbit polyclonal antibody for immunohistochemistry (paraffin) and ELISA. Targets human chorionic somatomammotropin hormone 1, a pregnancy-specific hormone involved in lactation, fetal growth, and metabolism. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	Placental lactogen (CSH1) Rabbit pAb
Remarks/Alias	PL; CSA; CS-1; CSMT; GHB3; hCS-1; hCS-A; Placental lactogen (CSH1)
Species	Human
Gene ID (Human)	1442
Gene ID	1442
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-217 of human Placental lactogen (CSH1) (NP_001308.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P0DML2
Protein Weight	25kDa
Shipping	Ice bag