

Anti-Plau Rabbit Polyclonal Antibody For Wb, Elisa - P00749

Item Number: CM0019504



Introduction

High-quality rabbit polyclonal antibody targeting human PLAU (urokinase-type plasminogen activator). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for protease and fibrinolysis research.

[Learn More](#)

Parameter	Value
Product Name	PLAU Rabbit pAb
Remarks/Alias	ATF; QPD; UPA; URK; u-PA; BDPLT5
Species	Human
Gene ID (Human)	5328
Gene ID	5328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-180 of human PLAU (NP_002649.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00749
Protein Weight	49kDa
Shipping	Ice bag