


Anti-Plaur Rabbit Polyclonal Antibody For Wb, Elisa - Q03405

Item Number: CM0013782



Introduction

High-quality anti-PLAUR rabbit polyclonal antibody for Western blot (WB) and ELISA, showing reactivity with human, mouse, and rat samples. Targets the urokinase receptor, critical for cell migration and tissue remodeling.

[Learn More](#)

Parameter	Value
Product Name	PLAUR Rabbit pAb
Remarks/Alias	CD87; UPAR; URKR; U-PAR; PLAUR
Species	Human
Gene ID (Human)	5329
Gene ID	5329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-305 of human PLAUR (NP_002650.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03405
Protein Weight	37kDa
Shipping	Ice bag