


Anti-Plc Beta 3 (Plcb3) Rabbit Monoclonal Antibody For Wb, Elisa - Q01970

Item Number: CM0002084



Introduction

High-specificity rabbit monoclonal antibody against human PLC beta 3 (PLCB3), suitable for Western blot and ELISA. Ideal for studying GPCR-mediated signaling, DAG/IP3 generation, and neutrophil function.

[Learn More](#)

Parameter	Value
Product Name	PLC beta 3 (PLCB3) Rabbit mAb
Remarks/Alias	SMDCD; PLC beta 3 (PLCB3)
Species	Human
GeneID (human)	5331
GeneID	5331
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01970
Protein Weight	139kDa
Shipping	Ice bag