


Anti-Plc Gamma 1 (Plcg1) Rabbit Monoclonal Antibody For Wb, Elisa - P19174

Item Number: CM0012266



Introduction

Rabbit monoclonal anti-PLCG1 antibody for WB and ELISA applications. Detects endogenous PLCγ1 in human, mouse, and rat samples. Suitable for signal transduction and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	PLC gamma 1 (PLCG1) Rabbit mAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1; PLC gamma 1 (PLCG1)
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 710-832 of human PLC gamma 1 (PLCG1) (P19174)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19174
Protein Weight	149 kDa
Shipping	Ice bag