


Anti-Plc Gamma 1 (Plcg1) Rabbit Polyclonal Antibody For Wb And Elisa - P19174

Item Number: CM0017177



Introduction

High-quality anti-PLC gamma 1 (PLCG1) rabbit polyclonal antibody validated for western blotting (WB) and ELISA. Recognizes human and mouse PLCG1 (149 kDa), a key phosphoinositide signaling enzyme involved in receptor tyrosine kinase pathways, actin reorganization, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	PLC gamma 1 (PLCG1) Rabbit pAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag