


Anti-Plc Gamma 2 (Plcg2) Polyclonal Antibody For Wb, If/Icc, Elisa - P16885

Item Number: CM0018744



Introduction

Rabbit polyclonal antibody targeting human PLC gamma 2 (PLCG2), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Crucial for transmembrane signaling and second messenger generation.

[Learn More](#)

Parameter	Value
Product Name	PLC gamma 2 (PLCG2) Rabbit pAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2; PLC gamma 2 (PLCG2)
Species	Human
GeneID (Human)	5336
GeneID	5336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PLC gamma 2 (PLCG2) (NP_002652.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag