


Anti-Plc Gamma 2 (Plcg2) Rabbit Monoclonal Antibody For Wb, Elisa - P16885

Item Number: CM0011487



Introduction

Rabbit monoclonal antibody targeting PLC gamma 2 (PLCG2), suitable for Western blot and ELISA applications. Cross-reacts with human and mouse species, enabling signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	PLC gamma 2 (PLCG2) Rabbit mAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2
Species	Human
Gene ID (Human)	5336
Gene ID	5336
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 701-800 of PLC gamma 2 (PLCG2) (NP_002652.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag