

Anti-Plcb1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nq66

Item Number: CM0022609

Introduction

Rabbit polyclonal antibody against PLCB1, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat PLCB1 with a predicted molecular weight of 139 kDa. Suitable for phospholipase C signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PLCB1 Rabbit pAb
Remarks/Alias	DEE12; PLC-I; EIEE12; PI-PLC; PLC154; PLCB1A; PLCB1B; PLC-154; PLC-beta-1; PLCB1
Species	Human
GeneID (human)	23236
GeneID	23236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 917-1216 of human PLCB1 (NP_056007.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ66
Protein Weight	139kDa
Shipping	Ice bag