

Anti-Plcb2 (Plc-Beta-2) Rabbit Polyclonal Antibody For Wb, Elisa - Q00722

Item Number: CM0029623



Introduction

Rabbit polyclonal antibody targeting PLCB2 (PLC-beta-2), validated for Western blot (WB) and ELISA. Detects a 134 kDa band in human, mouse, and rat samples. Suitable for studying GPCR signaling and phospholipase C pathways.

[Learn More](#)

Parameter	Value
Product Name	PLCB2 Rabbit pAb
Remarks/Alias	PLC-beta-2; PLCB2
Species	Human
GeneID (human)	5330
GeneID	5330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human PLCB2 (NP_004564.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00722
Protein Weight	134kDa
Shipping	Ice bag