

Anti-Plcd4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Brc7

Item Number: CM0025104

Introduction

PLCD4 Rabbit pAb detects human, mouse, and rat phospholipase C delta-4, suitable for Western blot and ELISA. Targets a phosphoinositide signaling enzyme generating DAG and IP3, involved in acrosome reaction and cell growth.

[Learn More](#)

Parameter	Value
Product Name	PLCD4 Rabbit pAb
Remarks/Alias	PLCD4
Species	Human
GeneID (Human)	84812
GeneID	84812
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PLCD4 (NP_116115.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRC7
Protein Weight	88kDa
Shipping	Ice bag