

Anti-Plch1 Polyclonal Antibody For Wb And Elisa - Q4Kwh8

Item Number: CM0018495

Introduction

PLCH1 Rabbit Polyclonal Antibody (Q4KWH8) validated for Western blot and ELISA applications. Specific for human PLCH1 with cross-reactivity to rat. Ideal for studying phospholipase C signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PLCH1 Rabbit pAb
Remarks/Alias	HPE14; PLCL3; PLC-L3; PLC eta 1; PLCH1
Species	Human
GeneID (Human)	23007
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1100-1400 of human PLCH1 (NP_055811.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q4KWH8
Protein Weight	189kDa
Shipping	Ice bag