

Anti-Plcxd2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q0Vaa5

Item Number: CM0017837

Introduction

Rabbit polyclonal antibody targeting human PLCXD2, validated for Western blot, immunofluorescence, and ELISA in human, mouse, and rat samples. Explore PI-PLC X domain-containing protein 2's role in lipid signaling.

[Learn More](#)

Parameter	Value
Product Name	PLCXD2 Rabbit pAb
Remarks/Alias	PLCXD2
Species	Human
GeneID (Human)	257068
GeneID	257068
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human PLCXD2 (NP_001172035.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q0VAA5
Protein Weight	35kDa
Shipping	Ice bag