


Anti-Pld2 Polyclonal Antibody For Wb, Elisa - O14939

Item Number: CM0013849



Introduction

Rabbit polyclonal antibody against human PLD2, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying phospholipase D2 in signal-induced cytoskeletal regulation and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	PLD2 Rabbit pAb
Remarks/Alias	PLD1C, PLD2
Species	Human
GeneID (Human)	5338
GeneID	5338
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human PLD2 (NP_002654.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14939
Protein Weight	106kDa
Shipping	Ice bag