

Anti-PLD3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iv08

Item Number: CM0023937



Introduction

Rabbit polyclonal antibody targeting human PLD3, validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Ideal for studying lysosomal exonuclease function, innate immunity, and nucleic-acid sensing.

[Learn More](#)

Parameter	Value
Product Name	PLD3 Rabbit pAb
Remarks/Alias	AD19, HUK4, HU-K4, SCA46, PLD3
Species	Human
GeneID (human)	23646
GeneID	23646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 279-490 of human PLD3 (NP_001026866)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q8IV08
Protein Weight	55kDa
Shipping	Ice bag