

Anti-PLD4 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Bz4

Item Number: CM0013659



Introduction

High-quality anti-PLD4 rabbit polyclonal antibody validated for WB and ELISA against human and rat PLD4. Ideal for nucleic-acid sensing and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	PLD4 Rabbit pAb
Remarks/Alias	C14orf175; PLD4
Species	Human
GeneID (Human)	122618
GeneID	122618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-200 of human PLD4 (NP_620145.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96BZ4
Protein Weight	56kDa
Shipping	Ice bag