


Anti-Pleckstrin Monoclonal Antibody For Wb, Elisa - P08567

Item Number: CM0003607



Introduction

High-quality rabbit monoclonal antibody specific for human Pleckstrin (p47), validated for Western blot and ELISA. Ideal for platelet signaling studies and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	Pleckstrin Rabbit mAb
Remarks/Alias	P47; PLEK1; Pleckstrin
Species	Human
GeneID (Human)	5341
GeneID	5341
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08567
Protein Weight	40 kDa
Shipping	Ice bag