


Anti-Plek Polyclonal Antibody For Wb, Elisa - P08567

Item Number: CM0029084



Introduction

Rabbit polyclonal antibody targeting human pleckstrin (PLEK/p47), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying protein kinase C substrate function in platelets and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	PLEK Rabbit pAb
Remarks/Alias	P47; PLEK1; PLEK
Species	Human
GeneID (human)	5341
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human PLEK (NP_002655.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08567
Protein Weight	40kDa
Shipping	Ice bag