

Anti-Plekhhb2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Cs7

Item Number: CM0023490



Introduction

Rabbit polyclonal antibody targeting human PLEKHB2 (evectin-2), suitable for WB, IHC-P, and ELISA. Detects a 25 kDa protein involved in retrograde endosome transport. Ideal for IVD research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	PLEKHB2 Rabbit pAb
Remarks/Alias	EVT2; PLEKHB2
Species	Human
GeneID (Human)	55041
GeneID	55041
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human PLEKHB2 (NP_001094093.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q96CS7
Protein Weight	25kDa
Shipping	Ice bag