

Anti-Plexin A1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9UIW2

Item Number: CM0011142



Introduction

Rabbit monoclonal antibody targeting human Plexin A1 (Q9UIW2). Suitable for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for semaphorin signaling, axon guidance, and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	Plexin A1 Rabbit mAb
Remarks/Alias	NOV; NOV; PLXN1; DWOPNED; PLEXIN-A1; Plexin A1
Species	Human
GeneID (Human)	5361
GeneID	5361
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human Plexin A1 (Q9UIW2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIW2
Protein Weight	211kDa
Shipping	Ice bag