

Anti-Plg Rabbit Polyclonal Antibody For Wb And Elisa - P00747

Item Number: CM0026408



Introduction

Rabbit polyclonal antibody against human plasminogen (PLG) for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for studying fibrinolysis, tissue remodeling, and coagulation pathways.

[Learn More](#)

Parameter	Value
Product Name	PLG Rabbit pAb
Remarks/Alias	HAE4; PLG
Species	Human
GeneID (Human)	5340
GeneID	5340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human PLG (NP_000292.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00747
Protein Weight	91kDa
Shipping	Ice bag