

Anti-Plin5 Polyclonal Antibody For Wb, Elisa - Q00G26

Item Number: CM0024217

Introduction

High-quality PLIN5 Rabbit pAb validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for lipid droplet and metabolism research. SWISS Q00G26.

[Learn More](#)

Parameter	Value
Product Name	PLIN5 Rabbit pAb
Remarks/Alias	MLDP; LSDA5; LSDP5; OXPAT
Species	Human
Gene ID (Human)	440503
Gene ID	440503
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 302-401 of human PLIN5 (NP_001013728.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q00G26
Protein Weight	51kDa
Shipping	Ice bag