

Anti-Plin5 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q8Bvz1

Item Number: CM0026704

Introduction

Anti-PLIN5 Rabbit Polyclonal Antibody validated for WB, IF-P, IHC-P, ELISA with Mouse and Rat reactivity. An essential tool for studying lipid droplet biology, fatty acid oxidation, and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	PLIN5 Rabbit pAb
Remarks/Alias	MLDP; Lsdp5; PAT-1; 2310076L09Rik
Species	Mouse
Gene ID (Human)	440503
Gene ID	66968
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 364-463 of mouse PLIN5 (NP_080150.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8BVZ1
Protein Weight	50 kDa/49 kDa
Shipping	Ice bag