


Anti-Plk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P53350

Item Number: CM0027898



Introduction

Rabbit polyclonal antibody targeting human PLK1, suitable for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat PLK1. Essential for studying mitotic progression, centrosome dynamics, and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	PLK1 Rabbit pAb
Remarks/Alias	PLK1, PLK, STPK13
Species	Human
GeneID (Human)	5347
GeneID	5347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 304-603 of human PLK1 (NP_005021.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53350
Protein Weight	68kDa
Shipping	Ice bag