


Anti-Plk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nyy3

Item Number: CM0025651



Introduction

Rabbit polyclonal antibody against human PLK2 (SNK/hPlk2), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying synaptic plasticity and centriole duplication.

[Learn More](#)

Parameter	Value
Product Name	PLK2 Rabbit pAb
Remarks/Alias	SNK; hSNK; hPlk2; PLK2
Species	Human
GeneID (Human)	10769
GeneID	10769
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 361-439 of human PLK2 (NP_006613.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYY3
Protein Weight	78kDa
Shipping	Ice bag