

Anti-Plk3 Polyclonal Antibody For Wb, Elisa - Q9H4B4

Item Number: CM0025665



Introduction

PLK3 Rabbit pAb validated for WB, ELISA in Human, Mouse, Rat. Recombinant fusion protein immunogen. Serine/threonine kinase PLK3 involved in cell cycle and stress response. Ideal for apoptosis and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	PLK3 Rabbit pAb
Remarks/Alias	CNK; FNK; PRK; PLK-3; PLK3
Species	Human
GeneID (Human)	1263
GeneID	1263
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 487-646 of human PLK3 (NP_004064.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4B4
Protein Weight	72kDa
Shipping	Ice bag