

## Anti-Plk4 Polyclonal Antibody For Wb And Elisa - O00444

Item Number: CM0025770



### Introduction

High-quality rabbit polyclonal antibody against PLK4 (Polo-like kinase 4) for Western blot (WB) and ELISA applications. Detects human, mouse, and rat PLK4 (~109 kDa).

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | PLK4 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SAK; STK18; MCCRP2; PLK4                                                                                          |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 10733                                                                                                             |
| Gene ID          | 10733                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 641-970 of human PLK4 (NP_055079.3) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O00444                                                                                                            |
| Protein Weight   | 109kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |