

Anti-PLP Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y342

Item Number: CM0014492

Introduction

High-quality rabbit polyclonal antibody against human PLLP, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for myelin membrane biogenesis and apical endocytosis studies.

[Learn More](#)

Parameter	Value
Product Name	PLLP Rabbit pAb
Remarks/Alias	PMLP, TM4SF11, PLLP
Species	Human
Gene ID (Human)	51090
Gene ID	51090
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PLLP (NP_057077.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y342
Protein Weight	20kDa
Shipping	Ice bag