

Anti-Plod1 Polyclonal Antibody For Wb, Elisa - Q02809

Item Number: CM0017300

Introduction

Rabbit polyclonal antibody targeting human PLOD1 (lysyl hydroxylase 1), suitable for Western blot and ELISA. Cross-reacts with mouse. Ideal for collagen biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PLOD1 Rabbit pAb
Remarks/Alias	LH; LH1; LLH; EDS6; PLOD; EDSKCL1; PLOD1
Species	Human
GeneID (human)	5351
GeneID	5351
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-727 of human PLOD1 (NP_000293.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q02809
Protein Weight	84kDa
Shipping	Ice bag