

Anti-Plod1 Rabbit Pab For Wb & Elisa - Q02809

Item Number: CM0030616

Introduction

Rabbit polyclonal antibody targeting human PLOD1 (lysyl hydroxylase 1) for Western blot and ELISA. Cross-reactive with mouse and rat. Suitable for collagen biosynthesis and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	PLOD1 Rabbit pAb
Remarks/Alias	LH; LH1; LLH; EDS6; PLOD; EDSKCL1; PLOD1
Species	Human
GeneID (Human)	5351
GeneID	5351
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-295 of human PLOD1 (NP_000293.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02809
Protein Weight	84kDa
Shipping	Ice bag