

Anti-Plod2/Lh2 Rabbit Monoclonal Antibody For Wb, Elisa - 000469

Item Number: CM0006977

Introduction

Recombinant rabbit monoclonal antibody targeting PLOD2/LH2 (lysyl hydroxylase 2), validated for WB and ELISA. Detects an 85kDa protein in human and rat samples. Essential for collagen cross-link and stability research.

[Learn More](#)

Parameter	Value
Product Name	PLOD2/LH2 Rabbit mAb
Remarks/Alias	LH2; TLH; BRKS2; PLOD2/LH2
Species	Human
GeneID (Human)	5352
GeneID	5352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-380 of human PLOD2/LH2 (NP_000926.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O00469
Protein Weight	85kDa
Shipping	Ice bag