

Anti-Plod2/Lh2 Rabbit Polyclonal Antibody For Wb And Elisa - 000469

Item Number: CM0027466

Introduction

PLOD2/LH2 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive across human, mouse, and rat. Targets procollagen-lysine,2-oxoglutarate 5-dioxygenase 2, a key enzyme in collagen cross-linking.

[Learn More](#)

Parameter	Value
Product Name	PLOD2/LH2 Rabbit pAb
Remarks/Alias	LH2; TLH; BRKS2; PLOD2/LH2
Species	Human
Gene ID (Human)	5352
Gene ID	5352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-380 of human PLOD2/LH2 (NP_000926.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	000469
Protein Weight	85kDa
Shipping	Ice bag