


Anti-Plod3/Lh3 Monoclonal Antibody For Wb, Elisa - O60568

Item Number: CM0002150



Introduction

Rabbit monoclonal antibody against human PLOD3 (LH3), validated for WB and ELISA. Recognizes human and mouse PLOD3, molecular weight 85kDa. Suitable for collagen modification research.

[Learn More](#)

Field	Value
Product Name	PLOD3 Rabbit mAb
Remarks/Alias	LH3; PLOD3
Species	Human
Gene ID (Human)	8985
Gene ID	8985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-306 of human PLOD3 (NP_001075.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60568
Protein Weight	85kDa
Shipping	Ice bag