

Anti-Plp1 Polyclonal Antibody For Wb, If-P, Elisa - P60201

Item Number: CM0013664



Introduction

Rabbit polyclonal antibody targeting human PLP1 (myelin proteolipid protein), suitable for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for myelin and neurological disorder research.

[Learn More](#)

Parameter	Value
Product Name	PLP1 Rabbit pAb
Remarks/Alias	PLP; PMD; HLD1; MMPL; SPG2; GPM6C; PLP/DM20; PLP1
Species	Human
GeneID	5354
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 102-148 of human PLP1 (NP_955772.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60201
Protein Weight	30kDa
Shipping	Ice bag