

Anti-Plp1 Polyclonal Antibody For Wb, Ihc, If, Elisa - P60201

Item Number: CM0023007



Introduction

Rabbit polyclonal antibody against human PLP1 (myelin proteolipid protein), validated for WB, IHC, IF, and ELISA. Reacts with human, mouse, and rat samples. Ideal for neuroscience and myelin research.

[Learn More](#)

Parameter	Description
Product Name	PLP1 Rabbit pAb
Remarks/Alias	PLP; PMD; HLD1; MMPL; SPG2; GPM6C; PLP/DM20; PLP1
Species	Human
GeneID	5354
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-150 of human PLP1 (NP_000524.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60201
Protein Weight	30kDa
Shipping	Ice bag