

Anti-Pls1 (Plastin-1) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14651

Item Number: CM0017634

Introduction

Rabbit polyclonal antibody against human plastin-1 (PLS1), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying actin bundling and stereocilia formation.

[Learn More](#)

Parameter	Value
Product Name	PLS1 Rabbit pAb
Remarks/Alias	DFNA76; PLS1
Species	Human
Gene ID (Human)	5357
Gene ID	5357
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human PLS1 (NP_001165783.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14651
Protein Weight	70kDa
Shipping	Ice bag