


Anti-Pls3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P13797

Item Number: CM0027378



Introduction

Anti-PLS3 Rabbit Polyclonal Antibody targets plastin-3 (T-plastin), an actin-bundling protein. Validated for WB, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for cytoskeleton and osteoporosis research.

[Learn More](#)

Parameter	Value
Product Name	PLS3 Rabbit pAb
Remarks/Alias	BMND18; T-plastin; PLS3
Species	Human
Gene ID (Human)	5358
Gene ID	5358
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P13797
Protein Weight	71kDa
Shipping	Ice bag