

Anti-Plscr3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nry6

Item Number: CM0015351

Introduction

Rabbit polyclonal antibody against human phospholipid scramblase 3 (PLSCR3). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Detects mitochondrial scramblase involved in cardiolipin translocation and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	PLSCR3 Rabbit pAb
Remarks/Alias	PLSCR3
Species	Human
GeneID	57048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human PLSCR3 (NP_065093.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRY6
Protein Weight	32kDa
Shipping	Ice bag