

## Anti-Plscr3 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Nry6

Item Number: CM0003241

### Introduction

High-specificity rabbit monoclonal antibody against human PLSCR3, validated for WB and ELISA. Cross-reacts with mouse and rat. Expected molecular weight ~32 kDa. Suitable for mitochondrial phospholipid scramblase research.

[Learn More](#)

| Parameter        | Value                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------|
| Product Name     | PLSCR3 Rabbit mAb                                                                               |
| Remarks/Alias    | PLSCR3                                                                                          |
| Species          | Human                                                                                           |
| Gene ID (Human)  | 57048                                                                                           |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PLSCR3 (Q9NRY6) |
| Source           | Rabbit                                                                                          |
| Category         | Monoclonal Antibodies                                                                           |
| Application      | WB, ELISA                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                               |
| SWISS            | Q9NRY6                                                                                          |
| Protein Weight   | 32kDa                                                                                           |
| Shipping         | Ice bag                                                                                         |