

Anti-Pltp Monoclonal Antibody For Wb, Ihc-P, Elisa - P55058

Item Number: CM0008407

Introduction

Recombinant rabbit anti-PLTP monoclonal antibody for WB, IHC-P, ELISA applications. Detects endogenous human, mouse, and rat phospholipid transfer protein (~55 kDa). Supports HDL remodeling and lipid transport research.

[Learn More](#)

Parameter	Value
Product Name	PLTP Rabbit mAb
Remarks/Alias	BPIFE; HDLCQ9; PLTP
Species	Human
Gene ID (Human)	5360
Gene ID	5360
Immunogen	Synthetic peptide corresponding to aa 200-300 of human PLTP (P55058)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P55058
Protein Weight	55kDa
Shipping	Ice bag