

Anti-Pltp Rabbit Polyclonal Antibody For Wb And Elisa - P55058

Item Number: CM0029391

Introduction

High-quality Anti-PLTP Rabbit Polyclonal Antibody validated for WB and ELISA applications, with cross-reactivity to human and mouse PLTP. Ideal for studying lipid transfer, HDL remodeling, and cholesterol metabolism.

[Learn More](#)

Parameter	Specification
Product Name	PLTP Rabbit pAb
Remarks/Alias	BPIFE; HDLCQ9; PLTP
Species	Human
GeneID (Human)	5360
GeneID	5360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 142-441 of human PLTP (NP_872617.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55058
Protein Weight	55kDa
Shipping	Ice bag