

Anti-Plvap Polyclonal Antibody For Wb, Elisa - Q9Bx97

Item Number: CM0018304

Introduction

Rabbit polyclonal antibody against human PLVAP (PV1), validated for WB and ELISA. Suitable for mouse and rat cross-reactivity. Ideal for studying endothelial fenestrae, vascular permeability, and diaphragms in microvascular biology.

[Learn More](#)

Parameter	Value
Product Name	PLVAP Rabbit pAb
Remarks/Alias	PV1; FELS; PV-1; gp68; DIAR10; PLVAP
Species	Human
Gene ID	83483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-442 of human PLVAP (NP_112600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BX97
Protein Weight	51kDa
Shipping	Ice bag