

Anti-Plxna3 Rabbit Polyclonal Antibody For Wb, Elisa - P51805

Item Number: CM0013182

Introduction

Anti-PLXNA3 rabbit polyclonal antibody validated for WB and ELISA applications. Detects human PLXNA3, a semaphorin receptor involved in axon guidance and neuronal development. Ideal for neuroscience and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PLXNA3 Rabbit pAb
Remarks/Alias	6.3; PLXN3; PLXN4; XAP-6; HSSEXGENE; PLXNA3
Species	Human
GeneID (Human)	55558
GeneID	55558
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 950-1150 of human PLXNA3 (NP_059984.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51805
Protein Weight	208kDa
Shipping	Ice bag